

## Nelfinavir Mesylate

Catalog No: #S6001


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Nelfinavir Mesylate
Brief Description	Inhibitors
Purification	98.00%
Target Name	HIV Protease inhibitor
Calculated MW	663.89
Formulation	C33H49N3O7S2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Cancer SMILES: CC1=C(C=CC=C1O)C(=O)N[C@@H](CSC2=CC=CC=C2)[C@@H](CN3C[C@H]4CCCC[C@H]4C[C@H]3C(=O)NC(C)(C)O)C(=O)O Pathways: Microbiology & virology; Proteases/Proteasome Receptor: HIV protease Boiling pt: Melting pt: Solubility: Ethanol: 100 mg/mL (150.62 mM); DMSO: 100 mg/mL (150.62 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Xiang T, et al. Y Lett. 2015 Aug 1;364(1):79-88.

Note: This product is for in vitro research use only