

## TAPI-1

Catalog No: #S6009


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	TAPI-1
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR inhibitor
Calculated MW	499.6
Formulation	C26H37N5O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other SMILES: CC(C)CC(CC(=O)NO)C(=O)N[C@@H](CC1=CC2=CC=CC=C2C=C1)C(=O)N[C@@H](C)C(=O)NCCN

[r,c:19,21,24,t:15,17] Pathways: Neuroscience Receptor: ADAM17/TACE Boiling pt: Melting pt: Solubility: DMSO: 99 mg/mL (198.15 mM) Ethanol: 99 mg/mL (198.15 mM) Water: 60 mg/mL warmed (120.09 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. M&uuml;llberg J, et al. J Immunol. 1995, 155(11), 5198-5205

Note: This product is for in vitro research use only