

Z-VAD-FMK

Catalog No: #S6013


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Z-VAD-FMK
Brief Description	Inhibitors
Purification	98.00%
Target Name	Caspase inhibitor
Calculated MW	467.49
Formulation	C22H30FN3O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:C[C@@H](C(=O)N[C@@H](CC(=O)OC)C(=O)CF)NC(=O)[C@H](C(C)C)NC(=O)OCC1=CC=CC=C1Pathways:Apoptosis; Proteases/ProteasomeReceptor:Pan-caspaseBoiling pt:Melting pt:Solubility:DMSO: 93 mg/mL (198.93 mM)Appearance:WhitePowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Slee EA, et al. Biochem J, 1996, 315 (Pt 1), 21-24.

Note: This product is for in vitro research use only