

Maraviroc

Catalog No: #S6016


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Maraviroc
Brief Description	Inhibitors
Purification	100.00%
Target Name	CCR antagonist
Calculated MW	513.67
Formulation	C29H41F2N5O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology
 SMILES: CC1=NN=C(N1C2C[C@H]3CC[C@@H](C2)N3CC[C@@H](C4=CC=CC=C4)NC(=O)C5CCC(CC5)(F)F)C(C)C
 Pathways: Immunology/Inflammation; Microbiology&virology
 Receptor: MIP-1α; MIP-1β; RANTES
 Boiling pt: Melting pt: Solubility: Ethanol: 100 mg/mL (194.67 mM); DMSO: 100 mg/mL (194.67 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Dorr P, et al. Antimicrob Agents ChemOther. 2005, 49(11), 4721-4732.

Note: This product is for in vitro research use only