

OSI-906 (Linsitinib)

Catalog No: #S6017


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	OSI-906 (Linsitinib)
Brief Description	Inhibitors
Purification	98.00%
Target Name	IGF-1R inhibitor; Others inhibitor
Calculated MW	421.49
Formulation	C26H23N5O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: CC1(O)CC(C1)C1=NC(=C2N1C=CN=C2N)C1=CC2=C(C=C1)C=CC(=N2)C1=CC=CC=C1
 Pathways: Others; Tyrosine Kinase/Adaptors Receptor: IGF-1R; Insulin Receptor; IRR
 Boling pt: Melting pt: Solubility: DMSO: 84 mg/mL (199.29 mM) Ethanol:

References

- Mulvihill MJ, et al. Future Med Chem, 2009, 1(6), 1153-1171.

Note: This product is for in vitro research use only