

MEK-4 Polyclonal Antibody

Catalog No: #41138

Package Size: #41138-1 50ul #41138-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MEK-4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	MEK-4 Polyclonal Antibody detects endogenous levels of MEK-4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human MEK-4 around the non-phosphorylation site of S80.
Target Name	MEK-4
Other Names	MAP2K4; JNKK1; MEK4; MKK4; PRKMK4; SEK1; SERK1; SKK1; Dual specificity mitogen-activated protein kinase kinase 4; MAP kinase kinase 4; MAPKK 4; JNK-activating kinase 1; MAPK/ERK kinase 4; MEK 4; SAPK/ERK kinase 1; SEK1; Stress-activated pro
Accession No.	Swiss-Prot: P45985NCBI Gene ID: 6416
Uniprot	P45985
GeneID	6416;
SDS-PAGE MW	44kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images

Western Blot analysis of 293 cells using MEK-4 Polyclonal Antibody



Note: This product is for in vitro research use only