

## MEL-1B-R Polyclonal Antibody

Catalog No: #41139

Package Size: #41139-1 50ul #41139-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MEL-1B-R Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	MEL-1B-R Polyclonal Antibody detects endogenous levels of MEL-1B-R protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human MEL-1B-R.
Target Name	MEL-1B-R
Other Names	MTNR1B; Melatonin receptor type 1B; Mel-1B-R; Mel1b receptor
Accession No.	Swiss-Prot: P49286NCBI Gene ID: 4544
Uniprot	P49286
GeneID	4544;
SDS-PAGE MW	40kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/5000.

Not yet tested in other applications.

## Images



Western Blot analysis of HepG2 cells using MEL-1B-R Polyclonal Antibody



Western Blot analysis of HELA A549 cells using MEL-1B-R Polyclonal Antibody

---

Note: This product is for in vitro research use only