

## Telotristat Etiprate (LX 1606 Hippurate)

Catalog No: #S6044


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Telotristat Etiprate (LX 1606 Hippurate)
Brief Description	Inhibitors
Purification	99.00%
Target Name	Hydroxylase activator
Calculated MW	754.15
Formulation	C27H26ClF3N6O3.C9H9NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism

systemSMILES:OC(=O)CNC(=O)C1=CC=CC=C1.CCOC(=O)[C@@H](N)CC1=CC=C(C=C1)C1=CC(O[C@@H](C2=CC=C(Cl)C=C2N2C=CC(C)=N2)C(F)(F)F)=NC(N)=N1 |r,c:9,11,23,25,38,42,45,51,54,t:7,21,28,33,35|Pathways:MetabolismReceptor:Tryptophan hydroxylaseBoiling pt:Melting pt:Solubility:DMSO:100 mg/mL (132.59 mM)Ethanol:

## References

1. Pappas C, et al. J Clin Oncol 27, 2009 (suppl; abstr e14555).

Note: This product is for in vitro research use only