

TAK-700 (Orteronel)


Catalog No: #S6050

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TAK-700 (Orteronel)
Brief Description	Inhibitors
Purification	/
Target Name	P450 inhibitor
Calculated MW	307.35
Formulation	C18H17N3O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: CNC(=O)C1=CC2=C(C=C1)C=C(C=C2)C3(CCN4C3=CN=C4)O Pathways: Metabolism Receptor: 17,20-lyase (Human) , 17,20-lyase (Rat) Boiling pt: Melting pt: Solubility: DMSO: 61 mg/mL (198.47 mM) Ethanol: 8 mg/mL (26.02 mM) Water:

Note: This product is for in vitro research use only