

## Turofexorate Isopropyl (XL335)

Catalog No: #S6053


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Turofexorate Isopropyl (XL335)
Brief Description	Inhibitors
Purification	99.33%
Target Name	FXR agonist
Calculated MW	438.47
Formulation	C25H24F2N2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular system/

Metabolism SMILES: CC(C)OC(=O)C1=CN(CC(C)(C)C2=CC=CC=C2)C(=O)C1=CC(F)=C(F)C=C1 Pathways: Metabolism Receptor: FXR

Boiling pt: Melting pt: Solubility: DMSO: 33 mg/mL (75.26 mM) Ethanol: 2 mg/mL (4.56 mM) Water:

## References

1. Flatt B, et al. J Med Chem, 2009, 52(4), 904-907.

Note: This product is for in vitro research use only