

Met Polyclonal Antibody

Catalog No: #41140

Package Size: #41140-1 50ul #41140-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Met Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Met Polyclonal Antibody detects endogenous levels of Met protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Met around the non-phosphorylation site of Y1234.
Target Name	Met
Other Names	MET; Hepatocyte growth factor receptor; HGF receptor; HGF/SF receptor; Proto-oncogene c-Met; Scatter factor receptor; SF receptor; Tyrosine-protein kinase Met
Accession No.	Swiss-Prot: P08581NCBI Gene ID: 4233
Uniprot	P08581
GeneID	4233;
SDS-PAGE MW	145kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using Met Polyclonal Antibody

Note: This product is for in vitro research use only