

Quisinostat (JNJ-26481585)

Catalog No: #S6055


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Quisinostat (JNJ-26481585)
Brief Description	Inhibitors
Purification	99.00%
Target Name	HDAC inhibitor
Calculated MW	394.476
Formulation	C21H26N6O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: CN1C=C(C2=CC=CC=C21)CNCC3CCN(CC3)C4=NC=C(C=N4)C(=O)NO
 Pathways: Chromatin/Epigenetic; DNA Damage/DNA Repair; NF-kB Receptor: HDAC1; HDAC10; HDAC11; HDAC2; HDAC3; HDAC4; HDAC5; HDAC8
 Boling pt: Melting pt: Solubility: DMSO: 79 mg/mL (200.26 mM) Ethanol:

References

1. Arts J, et al. Clin Cancer Res, 2009, 15(22), 6841-6851.

Note: This product is for in vitro research use only