

Thiamet G

Catalog No: #S6056


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thiamet G
Brief Description	Inhibitors
Purification	/
Target Name	Autophagy inhibitor
Calculated MW	248.3
Formulation	C9H16N2O4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:[H][C@]12O[C@H](CO)[C@@H](O)[C@H](O)[C@@]1([H])N=C(NCC)S2

[r,t:13|Pathways:AutophagyReceptor:O-GlcNAcaseBoling pt:Melting pt:Solubility:DMSO:50 mg/mL (201.36 mM)Ethanol:12 mg/mL (48.32 mM)Water50 mg/mL (201.36 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only