

Brefeldin A

Catalog No: #S6062


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Brefeldin A
Brief Description	Inhibitors
Purification	99.00%
Target Name	ATPase inhibitor
Calculated MW	280.36
Formulation	C16H24O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cardiovascular

SystemSMILES:C[C@H]1CCC/C=C/[C@@H]2C[C@@H](C[C@H]2[C@@H]/C=C/C(=O)O1)O

Pathways:Membrane transporter/Ion channelReceptor:ATPase (HCT 116)Boiling pt:Melting pt:Solubility:DMSO: 4 mg/mL (14.26 mM)Appearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Zhu JW, et al. Bioorg Med Chem. 2000, 8(2), 455-463.

Note: This product is for in vitro research use only