

## Mfn2 Polyclonal Antibody

Catalog No: #41141

Package Size: #41141-1 50ul #41141-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Mfn2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA IF
Species Reactivity	Hu Ms Rt
Specificity	Mfn2 Polyclonal Antibody detects endogenous levels of Mfn2 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human Mfn2.
Target Name	Mfn2
Other Names	MFN2; CPRP1; KIAA0214; Mitofusin-2; Transmembrane GTPase MFN2
Accession No.	Swiss-Prot: O95140NCBI Gene ID: 9927
Uniprot	O95140
GeneID	9927;
SDS-PAGE MW	86kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

IF1:50-1:200

Not yet tested in other applications.

## Images



Western Blot analysis of COLO205 cells using Mfn2 Polyclonal Antibody

---

Note: This product is for in vitro research use only