

BGT226 (NVP-BGT226)

Catalog No: #S6072


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	BGT226 (NVP-BGT226)
Brief Description	Inhibitors
Purification	98.00%
Target Name	PI3K inhibitor
Calculated MW	650.6
Formulation	C28H25F3N6O2.C4H4O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:OC(=O)C=C/C(O)=O.COC1=CC=C(C=N1)C1=CC=C2N=CC3=C(N(C(=O)N3C)C3=CC=C(N4CCNCC4)C(=C3)C(F)(F)F)C2=C1
 |c:11,13,20,41,49,t:9,16,18,22,30,32|Pathways:PI3K/Akt/mTOR signalingReceptor:PI3K α ; PI3K β ; PI3K γ Boling pt:Melting pt:Solubility:DMSO:30
 mg/mL (46.11 mM)Ethanol:

References

1. Markman B, et al. Ann Oncol, 2012, 23(9), 2399-2408.

Note: This product is for in vitro research use only