

Miliciclib (PHA-848125)

Catalog No: #S6081


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Miliciclib (PHA-848125)
Brief Description	Inhibitors
Purification	/
Target Name	CDK inhibitor; Trk receptor inhibitor
Calculated MW	460.57
Formulation	C25H32N8O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: CNC(=O)C1=NN(C)C2=C1C(C)(C)CC1=CN=C(NC3=CC=C(C=C3)N3CCN(C)CC3)N=C21

|c:8,22,24,t:4,15,17,20,35|Pathways: Cell Cycle/Checkpoint; Tyrosine Kinase/AdaptorsReceptor: CDK2/CyclinA; CDK4/CyclinD1; CDK5/p35;

CDK7/CyclinH; TrkA Boling pt: Melting pt: Solubility: DMSO: 92 mg/mL (199.75 mM) Ethanol:

Note: This product is for in vitro research use only