

Clopamide


Catalog No: #S6086

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Clopamide
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	345.84
Formulation	C ₁₄ H ₂₀ ClN ₃ O ₃ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: C[C@@H]1CCC[C@@H](N1NC(=O)C2=CC(=C(C=C2)Cl)S(=O)(=O)N)C
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO > 10mM
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

Note: This product is for in vitro research use only