

MGST2 Polyclonal Antibody

Catalog No: #41143

Package Size: #41143-1 50ul #41143-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MGST2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	MGST2 Polyclonal Antibody detects endogenous levels of MGST2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human MGST2.
Target Name	MGST2
Other Names	MGST2; GST2; Microsomal glutathione S-transferase 2; Microsomal GST-2; Microsomal GST-II
Accession No.	Swiss-Prot: Q99735NCBI Gene ID: 4258
Uniprot	Q99735
GeneID	4258;
SDS-PAGE MW	17kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 293T cells using MGST2 Polyclonal Antibody

Note: This product is for in vitro research use only