

## Cabotegravir (GSK744, GSK1265744)

Catalog No: #S6098


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cabotegravir (GSK744, GSK1265744)
Brief Description	Inhibitors
Purification	99.00%
Target Name	HIV Protease inhibitor
Calculated MW	405.35
Formulation	C19H17F2N3O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: C[C@H]1CO[C@H]2N1C(=O)C3=C(C(=O)C(=CN3C2)C(=O)NCC4=C(C=C(C=C4)F)F)O Pathways: Microbiology & virology; Proteases/Proteasome Receptor: HIV integrase Boiling pt: Melting pt: Solubility: DMSO: 38 mg/mL warmed (93.74 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Yoshinaga T, et al. Antimicrob Agents Chem Others. 2015, 59(1), 397-406.

Note: This product is for in vitro research use only