

Entospletinib (GS-9973)

Catalog No: #S6101

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Entospletinib (GS-9973)
Brief Description	Inhibitors
Purification	99.23%
Target Name	Syk inhibitor
Calculated MW	411.46
Formulation	C23H21N7O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:CancerSMILES:C1CN(CCO1)C1=CC=C(NC2=NC(=CN3C=CN=C23)C2=CC=C3C=NNC3=C2)C=C1

|c:14,17,27,31,34,t:7,9,12,19,23,25|Pathways:AngiogenesisReceptor:SykBoling pt:Melting pt:Solubility:DMSO: 43 mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only