

Darapladib (SB-480848)

Catalog No: #S6109


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Darapladib (SB-480848)
Brief Description	Inhibitors
Purification	99.00%
Target Name	Phospholipase inhibitor
Calculated MW	666.77
Formulation	C36H38F4N4O2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory system/ Cardiovascular system/

MetabolismSMILES:CCN(CC)CCN(CC1=CC=C(C=C1)C1=CC=C(C=C1)C(F)(F)F)C(=O)CN1C2=C(CCC2)C(=O)N=C1SCC1=CC=C(F)C=C1Pathway
s:MetabolismReceptor:Lp-PLA2Boiling pt:Melting pt:Solubility:DMSO:100 mg/mL warmed (149.97 mM)Ethanol:100 mg/mL (149.97 mM)Water:

References

1. Blackie JA, et al. Bioorg Med Chem Lett. 2003, 13(6), 1067-1070.

Note: This product is for in vitro research use only