

EX 527 (Selisistat)

Catalog No: #S6111

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	EX 527 (Selisistat)
Brief Description	Inhibitors
Purification	99.00%
Target Name	Sirtuin inhibitor
Calculated MW	248.71
Formulation	C13H13ClN2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: NC(=O)C1CCCC2=C1NC1=C2C=C(Cl)C=C1
Pathways: Chromatin/Epigenetic; DNA Damage/DNA Repair
Receptor: SIRT1
Boiling pt: Melting pt: Solubility: DMSO: 50 mg/mL (201.03 mM) Ethanol: 18 mg/mL (72.37 mM) Water:

References

1. Solomon JM, et al. Mol Cell Biol, 2006, 26(1), 28-38.

Note: This product is for in vitro research use only