

Skepinone-L

Catalog No: #S6130


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Skepinone-L
Brief Description	Inhibitors
Purification	98.00%
Target Name	p38 MAPK inhibitor
Calculated MW	425.42
Formulation	C24H21F2NO4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: OC[C@@H](O)COC1=CC=C2CCC3=CC(NC4=CC=C(F)C=C4F)=CC=C3C(=O)C2=C1

[r,c:21,24,26,32,t:6,8,12,16,18] Pathways: MAPK Receptor: p38α Boiling pt: Melting pt: Solubility: DMSO: 85 mg/mL (199.8 mM) Ethanol: 85 mg/mL (199.8 mM) Water

References

1. Koeberle SC, et al. Nat Chem Biol. 2011, 8(2), 141-143.

Note: This product is for in vitro research use only