

Pimasertib (AS-703026)

Catalog No: #S6131


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pimasertib (AS-703026)
Brief Description	Inhibitors
Purification	98.15%
Target Name	MEK inhibitor
Calculated MW	431.2
Formulation	C15H15FIN3O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: C1=CC(=C(C=C1)F)NC2=C(C=CN=C2)C(=O)NC[C@@H](CO)O
 Pathways: MAPK Receptor: MEK1/2 (MM cell line)
 Boiling pt: Melting pt: Solubility: DMSO: 86 mg/mL (199.44 mM)
 Appearance: Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kim K, et al. Br J Haematol, 2010, 149(4), 537-549.

Note: This product is for in vitro research use only