

## ENMD-2076 L-(+)-Tartaric acid

Catalog No: #S6132


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	ENMD-2076 L-(+)-Tartaric acid
Brief Description	Inhibitors
Purification	/
Target Name	Aurora Kinase inhibitor; c-RET inhibitor; FLT inhibitor; Src inhibitor; VEGFR inhibitor
Calculated MW	525.56
Formulation	C25H31N7O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cancer  
 SMILES: O[C@H]([C@@H](O)C(O)=O)C(O)=O.CN1CCN(CC1)C1=NC(\C=C\C2=CC=CC=C2)=NC(NC2=NNC(C)=C2)=C1  
 [r,c:24,26,28,36,38,t:17,22,32] Pathways: Angiogenesis Receptor: FLT3, RET, Aurora A, VEGFR3/FLT4, Src  
 Boling pt: Melting pt: Solubility: DMSO: 100 mg/mL warmed (190.27 mM) Ethanol:

Note: This product is for in vitro research use only