

Remodelin hydrobromide

Catalog No: #S6133


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Remodelin hydrobromide
Brief Description	Inhibitors
Purification	99.00%
Target Name	RAAS inhibitor
Calculated MW	363.28
Formulation	C15H15BrN4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer SMILES: Br.N#CC1=CC=C(C=C1)C1=CSC(NN=C2CCCC2)=N1 Pathways: Endocrinology/Hormones Receptor: acetyl-transferase

NAT10 Boling pt: Melting pt: Solubility: DMSO: 72 mg/mL warmed (198.19 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Larriou D, et al. Science. 2014, 344(6183), 527-532.

Note: This product is for in vitro research use only