

Icotinib

Catalog No: #S6153


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Icotinib
Brief Description	Inhibitors
Purification	98.00%
Target Name	EGFR inhibitor
Calculated MW	391.42
Formulation	C22H21N3O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:C#CC1=CC(NC2=NC=NC3=CC=C(OCCOCCOCCO4)C=C23)=CC=C1Pathways:Angiogenesis; JAK/STAT signaling; Tyrosine Kinase/AdaptorsReceptor:EGFRBoiling pt:Melting pt:Solubility:DMSO:78 mg/mL (199.27 mM)Ethanol:7 mg/mL (17.88 mM)Water:

References

1. Tan F, et al. Lung Y, 2012, 76(2), 177-182.

Note: This product is for in vitro research use only