

## Octanoic acidCPI-613

Catalog No: #S6157


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Octanoic acidCPI-613
Brief Description	Inhibitors
Purification	99.00%
Target Name	Dehydrogenase inhibitor
Calculated MW	388.59
Formulation	C22H28O2S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism systemSMILES:OC(=O)CCCC(CCSCC1=CC=CC=C1)SCC1=CC=CC=C1

[c:14,16,23,25,t:12]Pathways:MetabolismReceptor:α-ketoglutarate dehydrogenase; PDHBoiling pt:Melting pt:Solubility:DMSO:78 mg/mL (200.72 mM)Ethanol:78 mg/mL (200.72 mM)Water:

## References

1. Zachar Z, et al. J Mol Med (Berl). 2011, 89(11), 1137-1148.

Note: This product is for in vitro research use only