

Salirasib

Catalog No: #S6163


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Salirasib
Brief Description	Inhibitors
Purification	98.00%
Target Name	Ras inhibitor
Calculated MW	358.54
Formulation	C22H30O2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: CC(=CCC/C(=C/CC/C(=C/CSC1=CC=CC=C1C(=O)O)/C)/C
 Pathways: MAPK Receptor; PPMTase
 Boiling pt: Melting pt: Solubility: DMSO > 10mM
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Marciano D, et al. J Med Chem. 1995, 38(8), 1267-1272.

Note: This product is for in vitro research use only