

Ciproxifan

Catalog No: #S6164


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ciproxifan
Brief Description	Inhibitors
Purification	98.00%
Target Name	Histamine Receptor antagonist
Calculated MW	386.4
Formulation	C16H18N2O2.C4H4O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:C1CC1C(=O)C2=CC=C(C=C2)OCCCC3=CN=CN3Pathways:GPCR/G Protein; Immunology/Inflammation; NeuroscienceReceptor:H3 receptorBoiling pt:Melting pt:Solubility:DMSO:54 mg/mL (139.75 mM)Ethanol:54 mg/mL (139.75 mM)Water:1 mg/mL (2.58 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ligneau X, et al. J Pharmacol Exp Ther, 1998, 287(2), 658-666.

Note: This product is for in vitro research use only