

## Telatinib

Catalog No: #S6166


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Telatinib
Brief Description	Inhibitors
Purification	99.00%
Target Name	c-Kit inhibitor; PDGFR inhibitor; VEGFR inhibitor
Calculated MW	409.83
Formulation	C <sub>20</sub> H <sub>16</sub> ClN <sub>5</sub> O <sub>3</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cancer SMILES: CNC(=O)C1=NC=CC(COC2=NN=C(NC3=CC=C(Cl)C=C3)C3=C2OC=C3)=C1

|c:6,21,24,28,30,t:4,11,13,16,18|Pathways: Angiogenesis; Tyrosine Kinase/Adaptors Receptor: c-Kit; PDGFR $\alpha$ ; VEGFR2; VEGFR3

Boling pt: Melting pt: Solubility: DMSO: 82 mg/mL (200.08 mM) Ethanol: 1 mg/mL (2.44 mM) Water:

## References

1. Steeghs N, et al. Clin Cancer Res, 2008, 14(11), 3470-3476.

Note: This product is for in vitro research use only