

Moxalactam Disodium

Catalog No: #S6179


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Moxalactam Disodium
Brief Description	Inhibitors
Purification	/
Target Name	Antibacterial inhibitor
Calculated MW	564.44
Formulation	C20H18N6Na2O9S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: [H][C@]12OCC(CSC3=NN=NN3C)=C(N1C(=O)[C@]2(NC(=O)[C@H](C(=O)O[Na])C1=CC=C(O)C=C1)OC)C(=O)O[Na] |r,c:9,13,34,t:7,29,31|Pathways: Microbiology&virology Receptor: Antibacterial Boiling pt: Melting pt: Solubility: DMSO: 38mg/mL Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

Note: This product is for in vitro research use only