

Saxagliptin

Catalog No: #S6203


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Saxagliptin
Brief Description	Inhibitors
Purification	/
Target Name	DPP-4 inhibitor
Calculated MW	315.41
Formulation	C18H25N3O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: [H][C@@]12C[C@]1([H])N([C@@H](C2)C#N)C(=O)[C@@H](N)C12CC3CC(CC(O)(C3)C1)C2

|r|Pathways: Proteases/Proteasome Receptor: DPP-4 Boiling pt: Melting pt: Solubility: DMSO: 63 mg/mL (199.74 mM) Ethanol: 24 mg/mL (76.09

mM) Water: 63 mg/mL (199.74 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

Note: This product is for in vitro research use only