

Posaconazole

Catalog No: #S6211


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Posaconazole
Brief Description	Inhibitors
Purification	99.93%
Target Name	Others inhibitor
Calculated MW	700.78
Formulation	C37H42F2N8O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: OtherSMILES: CC[C@@H]([C@H](C)O)N1C(=O)N(C=N1)C2=CC=C(C=C2)N3CCN(CC3)C4=CC=C(C=C4)OC[C@H]5C[C@](OC5)(CN6C=NC=N6)C7=C(C=C(C=C7)F)F Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 100 mg/mL (142.69 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Benaim G, et al. J Med Chem. 2006, 49(3), 892-389.

Note: This product is for in vitro research use only