

Abiraterone Acetate

Catalog No: #S6215


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Abiraterone Acetate
Brief Description	Inhibitors
Purification	99.89%
Target Name	P450 inhibitor
Calculated MW	391.55
Formulation	C ₂₆ H ₃₃ NO ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

systemSMILES:CC(=O)O[C@H]1CC[C@@]2([C@H]3CC[C@]4([C@H]([C@@H]3CC=C2C1)CC=C4C5=CN=CC=C5)C)C

Pathways: Metabolism Receptor: CYP17

Boling pt: Melting pt: Solubility: Ethanol: 28 mg/mL (71.51 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Pinto-Bazurco Mendieta MA, et al. J Med Chem. 2008, 51(16), 5009-5018.

Note: This product is for in vitro research use only