

Abiraterone

Catalog No: #S6216


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Abiraterone
Brief Description	Inhibitors
Purification	99.87%
Target Name	P450 inhibitor
Calculated MW	349.51
Formulation	C24H31NO
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism

systemSMILES:C[C@]12CC[C@@H](CC1=CC[C@@H]3[C@@H]2CC[C@]4([C@H]3CC=C4C5=CN=CC=C5)C)O

Pathways:MetabolismReceptor:CY P17

Boling pt:Melting pt:Solubility:Ethanol: 0.2 mg/mL (0.57 mM); Water: 0.02 mg/mL (0.05 mM); DMSO: 0.1 mg/mL (0.28 mM)

Appearance:WhitePowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Attard G, et al. J Clin Oncol. 2008, 26(28), 4563-4571.

Note: This product is for in vitro research use only