

Irinotecan

Catalog No: #S6228


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Irinotecan
Brief Description	Inhibitors
Purification	99.00%
Target Name	Topoisomerase inhibitor
Calculated MW	586.68
Formulation	C33H38N4O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer
 SMILES: CCC1=C2C=C(C=CC2=NC3=C1CN4C3=CC5=C(C4=O)COC(=O)[C@@]5(CC)O)OC(=O)N6CCC(CC6)N7CCCCC7
 Pathways: DNA Damage/DNA Repair
 Receptor: Topo I
 Boiling pt: Melting pt: Solubility: DMSO: 7 mg/mL (11.93 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Pavillard V, et al. Y ChemOthers Pharmacol. 2002, 49(4), 329-335.

Note: This product is for in vitro research use only