

## Tebipenem Pivoxil

Catalog No: #S6231


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Tebipenem Pivoxil
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	497.63
Formulation	C22H31N3O6S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:[H][C@]12[C@@H](C)C(SC3CN(C3)C3=NCCS3)=C(N1C(=O)[C@]2([H])[C@@H](C)O)C(=O)OCOC(=O)C(C)(C)C  
 |r,c:16,t:11|Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:99 mg/mL (198.94 mM)Ethanol:87 mg/mL (174.82 mM)Water

## References

1. Kato K, et al. Mol Pharm, 2010, 7(5), 1747-1756.

Note: This product is for in vitro research use only