

Lapatinib (GW-572016) Ditosylate

Catalog No: #S6235


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lapatinib (GW-572016) Ditosylate
Brief Description	Inhibitors
Purification	99.00%
Target Name	EGFR inhibitor; Raf inhibitor; Src inhibitor
Calculated MW	925.46
Formulation	C29H26ClFN4O4S2C7H8O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer
 SMILES: CC1=CC=C(C=C1)S(O)(=O)O.CC1=CC=C(C=C1)S(O)(=O)O.CS(=O)(=O)CCNCC1=CC=C(O1)C1=CC2=C(C=C1)N=CN=C2NC1=CC(Cl)=C(OCC2=CC(F)=CC=C2)C=C1
 Pathways: Angiogenesis; MAPK; JAK/STAT signaling; Tyrosine Kinase/Adaptors Receptor: EGFR; HER2/ErbB2; HER2/ErbB2; C-Raf-1; c-Src
 Boiling pt: Melting pt: Solubility: DMSO: 100 mg/mL (108.05 mM)
 Appearance: Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Rusnak DW, et al. Mol Y Ther, 2001, 1(2), 85-94.

Note: This product is for in vitro research use only