

Ozagrel

Catalog No: #S6236


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ozagrel
Brief Description	Inhibitors
Purification	99.00%
Target Name	Thrombin inhibitor
Calculated MW	228.25
Formulation	C13H12N2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:C1=CC(=CC=C1CN2C=CN=C2)/C=C/C(=O)OPathways:OthersReceptor:Thromboxane A2 synthetaseBoiling pt:Melting pt:Solubility:DMSO: 46 mg/mL (201.53 mM)Appearance:WhitePowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ishitsuka Y, et al. J Pharmacol Sci, 2009, 111(2), 211-215.

Note: This product is for in vitro research use only