

Trelagliptin

Catalog No: #S6237


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trelagliptin
Brief Description	Inhibitors
Purification	98.00%
Target Name	DPP-4 inhibitor
Calculated MW	357.38
Formulation	C18H20FN5O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:CN1C(=O)C=C(N(C1=O)CC=C(C=CC(=C2)F)C#N)N3CCC[C@H](C3)NPathways:Proteases/ProteasomeReceptor:DPP-4Bolin g pt:Melting pt:Solubility:Ethanol: 3 mg/mL warmed (8.39 mM); DMSO: 71 mg/mL (198.66 mM)Appearance:WhitePowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Inagaki N, et al. Lancet Diabetes Endocrinol. 2014, 2(2), 125-132.

Note: This product is for in vitro research use only