

BAPTA-AM

Catalog No: #S6245


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	BAPTA-AM
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	764.68
Formulation	C34H40N2O18
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:CC(=O)OCOC(=O)CN(CC(=O)OCOC(C)=O)C1=CC=CC=C1OCCOC1=CC=CC=C1N(CC(=O)OCOC(C)=O)CC(=O)OCOC(C)=OP
 athways:OthersReceptor:OthersBoling pt:Melting pt:Solubility:DMSO: 20 mg/mL (26.15 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wie MB, et al. Prog Neuropsychopharmacol Biol Psychiatry. 2001, 25(8), 1641-1659.

Note: This product is for in vitro research use only