

Cobicistat (GS-9350)

Catalog No: #S6246


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cobicistat (GS-9350)
Brief Description	Inhibitors
Purification	98.00%
Target Name	P450 inhibitor
Calculated MW	776.02
Formulation	C40H53N7O5S2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

systemSMILES:CC(C)C1=NC(CN(C)C(=O)N[C@@H](CCN2CCOCC2)C(=O)N[C@H](CC[C@H](CC2=CC=CC=C2)NC(=O)OCC2=CN=CS2)CC2=CC=CC=C2)=CS1 |r,c:32,34,44,51,53,55,t:3,30,42,49|Pathways:MetabolismReceptor:CYP3ABolig pt:Melting pt:Solubility:DMSO:100 mg/mL (128.86 mM)Ethanol:100 mg/mL (128.86 mM)Water

References

- Xu L, et al. ACS Med. Chem. Lett, 2010, 1(5), 209-213.

Note: This product is for in vitro research use only