


## Description

Product Name	XL888
Brief Description	Inhibitors
Purification	98%
Target Name	HSP inhibitor
Calculated MW	503.64
Formulation	C29H37N5O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:cancerSMILES:CC[C@@H](C)NC1=C(C=C(C)C(=C1)C(=O)NC1CC2CCC(C1)N2C1=NC=C(C=C1)C(=O)C1CC1)C(N)=O

[r,c:10,28,30,t:5,7,26]Pathways:Cytoskeletal Signaling; MetabolismReceptor:HSP90Boiling pt:Melting pt:Solubility:10 mM in

DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only