

Ipatasertib (GDC-0068)

Catalog No: #S6252


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ipatasertib (GDC-0068)
Brief Description	Inhibitors
Purification	99.00%
Target Name	Akt inhibitor
Calculated MW	458
Formulation	C24H32ClN5O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: CC(C)NC[C@@H](C(=O)N1CCN(CC1)C1=C2[C@H](C)C[C@@H](O)C2=NC=N1)C1=CC=C(Cl)C=C1

|r,c:15,23,25,33,t:28,30|Pathways: Cytoskeletal Signaling; PI3K/Akt/mTOR signaling Receptor: Akt1; Akt2; Akt3 Boling pt: Melting pt: Solubility: DMSO: 92 mg/mL (200.87 mM) Ethanol: 92 mg/mL (200.87 mM) Water

References

1. Lin K. Cancer Res, 2011, 71(8 Supplement), abstract DDT02-01.

Note: This product is for in vitro research use only