

## Bromosporine

Catalog No: #S6255


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Bromosporine
Brief Description	Inhibitors
Purification	99.36%
Target Name	CDK inhibitor; Epigenetic Reader Domain inhibitor
Calculated MW	404.44
Formulation	C17H20N6O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:othersSMILES:CCOC(=O)NC1=CC(=NN2C(C)=NN=C12)C1=CC(NS(C)(=O)=O)=C(C)C=C1 |c:8,12,28,t:6,14,18,25|Pathways:Cell Cycle/Checkpoint; Chromatin/EpigeneticReceptor:CECR2; BRD2; BRD4; BRD9Boiling pt:Melting pt:Solubility:DMSO: 51.7

mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only