


## Description

|                   |   |
|-------------------|---|
| Product Name      | YYA-021   |
| Brief Description | Inhibitors                                      |
| Purification      | 98.00%  |
| Target Name       | gp120/CD4 inhibitor                             |
| Calculated MW     | 317.43  |
| Formulation       | C18H27N3O2                                      |
| Storage           | 3 years -20°C powder; 2 years -80°C in solvent; |

## Images



## Product Description

Research Area: Other  
 SMILES: CC1=CC=C(C=C1)NC(=O)C(=O)NC2CC(NC(C2)(C)C)(C)C  
 Pathways: Immunology/Inflammation  
 Receptor: Gp120  
 Boiling pt: Melting pt: Solubility: DMSO: 0.01 mg/mL (0.03 mM) Ethanol:

## References

- Otsuki H, et al. J Gen Virol. 2013, 94(Pt 12), 2710-2716.

Note: This product is for in vitro research use only