

## MNDA Polyclonal Antibody

Catalog No: #41159

Package Size: #41159-1 50ul #41159-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	MNDA Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA IF
Species Reactivity	Hu Rt Ms
Specificity	MNDA Polyclonal Antibody detects endogenous levels of MNDA protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human MNDA.
Target Name	MNDA
Other Names	MNDA; Myeloid cell nuclear differentiation antigen
Accession No.	Swiss-Prot: P41218NCBI Gene ID: 4332
Uniprot	P41218
GeneID	4332;
SDS-PAGE MW	46kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

## Images



Western Blot analysis of 3T3 cells using MNDA Polyclonal Antibody

---

Note: This product is for in vitro research use only