

Bosentan

Catalog No: #S6264


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bosentan
Brief Description	Inhibitors
Purification	98.00%
Target Name	Endothelin Receptor antagonist
Calculated MW	551.61
Formulation	C27H29N5O6S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: COC1=CC=CC=C1OC1=C(NS(=O)(=O)C2=CC=C(C=C2)C(C)(C)C)N=C(N=C1OCCO)C1=NC=CC=N1
 |c:4,6,10,18,20,27,29,38,40,t:2,16,36|Pathways:GPCR/G ProteinReceptor:ET-A; ET-BBoiling pt:Melting pt:Solubility:DMSO:100 mg/mL (181.28 mM)Ethanol:3 mg/mL (5.43 mM)Water

References

1. Serafim KG, et al. Naunyn Schmiedebergs Arch Pharmacol. 2015 Nov;388(11):1211-21.

Note: This product is for in vitro research use only